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EN User manual for Alginate dressing

Description

Alginate dressing is a primary wound dressing made from the calcium salt of alginate acid. It is prepared as a textile fibre held together by a special needle bonding process, and presented as a “rope” for packing for cavities, and a flat non-woven pad for application to larger open wounds.

Indications

Alginate dressing may be applied to moderately to heavily exuding wounds such as partial-thickness burns, donor sites, leg ulcers, pressure ulcers, cavity wounds and most other granulating wounds. It achieves haemostasis.

Contra-indications

We recommend not using alginate dressing on dry wounds or wounds covered with hard, black necrotic tissue.

Application method

Alginate dressing, chosen to be slightly larger than the area of the wound, is placed in intimate contact with the wound base. Any overlap may be cut away or folded back over the wound surface before the application of a sterile secondary dressing held in place with surgical tape or a bandage as needed. The secondary dressing should be appropriate to the condition of the wound. If large quantities of exudate are anticipated, a simple absorbent dressing pad may be required; but as the wound heals and less exudate is produced, a thinner, low-tack pad with a plastic film may help to conserve moisture and prevent the wound drying out too quickly. Deeper cavity wounds or sinuses may be dressed with alginate rope, which should be placed gently in position but not packed in too tightly.

Frequency of dressing changes

The interval between dressing changes will depend entirely upon the state of the wound and the nature of the secondary dressing. The interval between changes may be extended to up to 7 days in some circumstances.

Warnings

Wounds that show signs of clinical infection may be dressed with this product but the dressing should be changed daily and the use of systemic antibiotic therapy considered.

Storage

Alginate dressings are stored in a well-ventilated clean room with relative humidity less than 80%, and no corrosive gases.

Other

STERILE Alginate dressings have been sterilized using Gamma radiation.


 Single sterile barrier system.

 Do not reuse product.

 Do not use if package has been compromised.

 Protect against direct sunlight.

 Keep dry.

 Store at a normal room temperature 5 - 35°C.

 Instructions for use.

DECIFERA is a registered trademark of Van Heek Medical. All rights reserved.

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