

BACTERIOSTATIC BLANKET With antibacterial and silver ion treatment



Perfect hygienic, decontaminable, heat insulated and highly comfortable blanket. Particularly useful for first aid, operating blocks and hemodialysis services.

Covers and protects the patient from cold. The bacteriostatic blanket keeps the patient warm thanks to the thermic performance of the HEATUP technical fibre. This whirl fibre regulates and restores the body heat. Very high mechanic and time resistance.

Bacteriostatic envelope made out of light and comfortable polyester coated with polyurethane. Fire-proof M1.

Ventilated : eases the air flow, thermic-variable coat, and adapts itself to the human warmth.

Waterproof : liquids and secretions proofed.

Bacteriostatic : Antibacterial and silver AG2+ ion treatment.

Finishing work : ultrasound welded

Dimensions : COUV1 : 110 x 190 cm or couv 2 : 140 x 200 cm

Color : green medical, navy blue, dark grey, dark blue.

Maintenance : with decontaminating spray or wipes.

We recommend to use the high surface bactericidal products of the ANIOS brand.

Washing 90°C – Drying 60°C- Test the drying before using, it took between 45 mn and 75 mn and depends of the machines.



Sarl GetC CLINIBED - 1, rue Leuridan-Noclain. 59150 Wattrelos – FRANCE

Tél : +33(0) 320-02-23-60 Fax : +33(0) 320-02-23-50

clinibed@clinibed.com - www.clinibed.com